

## **Highlights of the Hitachi AIRIS Mate .2T Open Bore MRI**

Considered #1 choice in the world for musculoskeletal and sports medicine professionals

Compact enough to be installed in rooms smaller than 25m squared

Includes a floating table that make positioning patients very easy, even when using a single finger.

### **Features**

- Whole-body, low-field MR Imaging
- Patient-friendly design - easy to maneuver
- Variety of fat-saturation techniques
- Low noise emissions - very quiet
- Considered the most cost-effective Open-Bore MRI

### **Specifications**

- Tesla Strength: 0.2T
- Dual-pillar design
- 210 degrees frontal opening
- 70 degrees rear opening
- Gradients: 15mTm; 30T/m/s
- Choice of patient table orientation, suites rural areas best
- Patient Table: Floating

### **Packages**

- Angiography (PC-angio, TOF-angio, TOF CE-MRA)
- Selective coil packages
- Fat-Saturation (DIXON/ STIR)
- Neuro Imaging
- Multiple Array

### **Coils**

- Bilateral Breast Coil
- Joint Coil S/M/L
- Multiple Array Head Coil
- Multiple Array Neck Coil
- QD Body Coil S/M/L
- QD Head Coil
- QD Knee Coil
- QD Wrist Coil
- Shoulder Coil

- TMJ Coils
- Vascular Head/ Neck Coil